


# ADJUSTABLE LATERAL POSITIONING SYSTEM




**CUSTOM  
SIZED  
SYSTEMS  
AVAILABLE**

- Center area is hole-free to eliminate ghosting
- 8 rods (2 sizes) for versatility in positioning
- Includes 4 acrylic blocks for finite positioning
- Radiolucent gel bumps included to simulate a breaking table
- Surgical bar rails available in either composite or aluminum material
- Standard ALPS comes configured for 22" or 24" wide tables



▲ Shown above with Travel Table, high speed drill, Bovie, and Flex-Light



Aluminum Rail Option ▲

▲ Composite Rail for 3D Navigation

Need to position a cadaveric specimen for a lateral procedure and a breaking OR table is not available? The ALPS by MEDSource is your solution. The ALPS allows you to simulate the proper patient positioning in a bioskills education environment.

MEDSource is your one-stop bioskills rental source. A variety of rental equipment (i.e. suction, towers, instruments, Bovie, lightsource, etc) is available.



Positioning tape and/or co-ban wrap are available for purchase. Don't forget to order PPE (gowns, masks, etc.) and Station kits (drapes, chux, irrigation bulb syringes, etc.) for your lab.

**TO LEARN MORE, CONTACT US TODAY!**

**MEDSource, Inc.**

7002 S Revere Parkway, Suite 25 | Centennial, Colorado 80112  
Ph: (303) 750-5357 | Fax: (303) 695-5276  
E-mail: [info@medsourcerental.com](mailto:info@medsourcerental.com)